

(NU3-6) Fundamental law of radioactive decay

Aim of experiment

Determination of the characteristic decay law for a short-lived radioactive source and its half life time.

Apparatus

GM tube counting station consists of GM counter – radioactive source- source holder- stop watch, source cabinet made of thick lead.

Theory of experiment

Each nucleus has a fixed probability of decaying per unit time. Nothing affects this probability (e.g., temperature, pressure, bonding environment, etc.) [Exception: very high pressure promotes electron capture slightly]

This is equivalent to saying that averaged over a large enough number of atoms the number of decays per unit time is proportional to the number of atoms present. Therefore in a closed system:

$$\frac{dN}{dt} \propto N(t)$$
$$\frac{dN}{dt} = -\lambda N(t) = A(t)$$

where A is the activity of the source after elapsed time, t .

and N = number of parent nuclei at time t and it is proportional to the activity A .

λ = decay constant = probability of decay per unit time. λ may have the units time^{-1} (e.g. s^{-1} , d^{-1} , m^{-1} , y^{-1})

To get time history of number of parent nuclei, integration of the above equation yields:

$$N(t) = N_0 e^{-\lambda t}$$

N_0 = initial number of parent nuclei at time $t = 0$.

At $t = T_{1/2}$, the half-life; the time it takes for half the nuclei to decay - i.e. for activity, A_0 , to decrease to $A_0/2$ is given by

$$T_{1/2} = \frac{\ln 2}{\lambda}$$




Figure 1 represents the relation between $A = \lambda N$ and time of decay.

At $t=0$, $N_0 = A_0/\lambda$, and at $t=T_{1/2}$, $N(T_{1/2}) = A_0/2\lambda$. One can obtain the half-life time and the decay constant of the decaying nucleus. If one plots the relation between $\ln(N)$ and the time of decay, t , a straight line with negative slope is obtained, from which λ is obtained.

$$\ln(N) = \ln(N_0) - \lambda t$$

Procedure

1. Set up the electronics as shown in Figure 2.
2. Place a radioactive source on a shelf of the GM stand at a suitable distance.
3. Switch on the power supply and let it warm for few minutes, then set the applied voltage at the tube operating voltage, V_{op}

4. Set the timer on 1 min and measure the background reading, N_{bg} c/min.
5. Start count, then start stop watch, $t=0$ min, and record the corresponding number of counts, N c/min.
6. For successive elapsed times, suitable for the source lifetime, record number of counts per minute in a table.
7. Draw the relation between $\ln(N \text{ c/min})$ and the elapsed time, t .
8. Draw the best line and find the slope, which equals the source disintegration constant, λ .

Results

$$N_{bg} = \text{c/min}$$

Time, t	N c/min	$N = N - N_{bg}$ c/min	$\ln(N \pm \sqrt{N})$

$$\lambda =$$



Figure 2. A schematic diagram of Geiger counting system