

(HT1-12) Boyle's Law of Gases

Aim of experiment

Determination of the atmospheric pressure.

Apparatus

Boyle's law apparatus – ruler – two stands.

Theory of experiment

Boyle's law states that for a given mass of gas maintained at constant temperature, the volume (V) is inversely proportional to the pressure (P), or $PV = \text{constant}$. A plot of $1/V$ against P will thus yield a straight line graph passing through the origin, *figure 1*.




Figure 1. The plot of $1/V$ against P

Now if H is the atmospheric pressure in cm of mercury, and h the difference in the mercury levels in tubes T_2 and T_1 , then the pressure in cm mercury is $P = H \pm h$. thus the plot of $1/V$ against h when extrapolated to cut the h - axis will locate the common zero, and the intercept, oA , is evidently the atmospheric pressure H , *figure 2*.





Figure 1. A schematic diagram of Boyle's law apparatus

Procedures

1. Start from the point at which the mercury level in the open end tube, T_2 , equal that of the closed end tube, T_1 and measure air volume contained in the graduated tube A, *figure 1*.

3. Move up the open end tube T_2 each 1 cm from the mercury level in the tube T_1 , which represents the height difference $h = C - B$ and record the corresponding volume.

Repeat the above steps two more times.

4. Record the results in table and find the average volume.

5. Draw a graph between $1/V_{av}$ on y-axis and h on the x-axis, a straight line is obtained, its intercept with $-ve$ x-axis represents the atmospheric pressure H .

Results

h (cm Hg)	V_1 (cm ³)	V_2 (cm ³)	V_3 (cm ³)	V_{av} (cm ³)	$1/V_{av}$ (cm ⁻³)

Density of mercury (ρ) = 1.335×10^3 Kg m⁻³

Acceleration of gravity (g) = 9.81 N Kg⁻¹

Atmospheric pressure in Pascal = $\rho g H =$